


Product Datasheet

Ly6g Antibody (Biotin) (orb1251121)

Description	Ly6g Antibody (Biotin)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC, FLISA
Target	Ly6g
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Ly-6G/Ly-6C Antibody [RB6-8C5] for use in flow cytometry, immunohistochemistry, and immunoprecipitation assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P35461
NCBI	P35461
Dilution Range	Ly-6G/Ly-6C Antibody [RB6-8C5] for use in flow cytometry, immunohistochemistry, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



BALB/c mouse bone marrow cells were stai...