


Product Datasheet

Kit Antibody (Biotin) (orb1251143)

Description	Kit Antibody (Biotin)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC, ICC, IHC-Fr, IP
Immunogen	Mouse IL-3-dependent bone marrow mast cells
Target	Kit
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD117 Antibody [2B8] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and immunoprecipitation assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P05532
NCBI	P05532
Dilution Range	CD117 Antibody [2B8] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



C57BL/6 mouse bone marrow cells were sta...