

Product Datasheet

LPAM-1 Antibody [DATK32] (orb1251146)

Description	LPAM-1 Antibody [DATK32]
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	AKR/Cum mouse T lymphoma cell line TK-1
Target	Itga4; Itgb7
Preservatives	Supplied in PBS
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at -20°C
Note	For research use only
Application notes	LPAM-1 Antibody [DATK32] for use in flow cytometry, immunohistochemistry, immunoprecipitation, and blocking assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P26011 , Q00651
NCBI	Q00651 , P26011
Dilution Range	LPAM-1 Antibody [DATK32] for use in flow cytometry, immunohistochemistry, immunoprecipitation, and blocking assays.
Expiration Date	12 months from date of receipt.