

Product Datasheet

MTRNR2L1 antibody (orb125123)

Description	Goat polyclonal antibody to MTRNR2L1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Humanin 1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Not yet tested - our routinely used western blotting protocol does not allow detection of proteins as small as the calculated size of 2.7kDa according to NP_001177381.1. Therefore we cannot recommend an optimal concentration and the antibody is an aspirin
Protein Sequence	TSEIDLPKRRRT
Clonality	Polyclonal
MW	2.69
Entrez	100462977
Dilution Range	ELISA: 1:16000
Expiration Date	12 months from date of receipt.