


Product Datasheet

Cd69 Antibody (FITC) (orb1251256)

| | |
|----------------------------|---|
| Description | Cd69 Antibody (FITC) |
| Species/Host | Hamster |
| Reactivity | Mouse |
| Conjugation | FITC |
| Tested Applications | FC, IHC-Fr, IP |
| Immunogen | Mouse dendritic epidermal cell line Y245 |
| Target | Cd69 |
| Preservatives | Supplied in PBS/NaN3 |
| Form/Appearance | Liquid |
| Concentration | 0.5 mg/mL |
| Storage | Store vial at 2-8°C |
| Note | For research use only |
| Application notes | CD69 Antibody [H1.2F3] for use in flow cytometry and immunohistochemistry assays. |
| Isotype | IgG |
| Clonality | Monoclonal |
| Uniprot ID | P37217 |
| NCBI | P37217 |
| Dilution Range | CD69 Antibody [H1.2F3] for use in flow cytometry and immunohistochemistry assays. |
| Expiration Date | 12 months from date of receipt. |



Con-A
stimulated
BALB/c
mouse
splenocyte...