

Product Datasheet

Itga5 Antibody (FITC) (orb1251267)

Description	Itga5 Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	C57BL/6 x A/J F1 mouse mast cell line MC/9
Target	Itga5
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD49e Antibody [5H10.27] for use in flow cytometry and immunohistochemistry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P11688
NCBI	P11688
Dilution Range	CD49e Antibody [5H10.27] for use in flow cytometry and immunohistochemistry assays.
Expiration Date	12 months from date of receipt.