


Product Datasheet

Cd40Ig Antibody (Biotin) (orb1251284)

Description	Cd40Ig Antibody (Biotin)
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC, ICC, IHC-Fr, IHC-P, IP, WB
Immunogen	Activated mouse Th1 clone D1.6
Target	Cd40Ig
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD154 Antibody [MR1] for use in flow cytometry, ELISA, immunohistochemistry, and immunoprecipitation assays.
Isotype	IgG
Clonality	Monoclonal
Uniprot ID	P27548
NCBI	P27548
Dilution Range	CD154 Antibody [MR1] for use in flow cytometry, ELISA, immunohistochemistry, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



PMA and Ionomycin-stimulated BALB/c mous...