


Product Datasheet

CD16/32 Antibody [93] (Biotin) (orb1251300)

Description	CD16/32 Antibody [93] (Biotin)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC
Immunogen	VH81X Tg B cells and T220 stromal cell line
Target	Fcgr2b; Fcgr3
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD16/32 Antibody [93] for use in flow cytometry, immunohistochemistry, and immunoprecipitation assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P08101 , P08508
NCBI	P08101 , P08508
Dilution Range	CD16/32 Antibody [93] for use in flow cytometry, immunohistochemistry, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



BALB/c mouse splenocytes were stained with...