

## Product Datasheet

### CD16/32 Antibody [93] (FITC) (orb1251301)

<b>Description</b>	CD16/32 Antibody [93] (FITC)
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	FC
<b>Immunogen</b>	VH81X Tg B cells and T220 stromal cell line
<b>Target</b>	Fcgr2b; Fcgr3
<b>Preservatives</b>	Supplied in PBS/NaN <sub>3</sub>
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.5 mg/mL
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	CD16/32 Antibody [93] for use in flow cytometry and immunohistochemistry assays.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P08101</a> , <a href="#">P08508</a>
<b>NCBI</b>	<a href="#">P08101</a> , <a href="#">P08508</a>
<b>Dilution Range</b>	CD16/32 Antibody [93] for use in flow cytometry and immunohistochemistry assays.
<b>Expiration Date</b>	12 months from date of receipt.