


Product Datasheet

Il2ra Antibody (Biotin) (orb1251320)

Description	Il2ra Antibody (Biotin)
Species/Host	Rat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	FC, ICC, IHC-Fr, IHC-P, WB
Immunogen	IL-2 dependent BALB/c helper T cell clone HT2
Target	Il2ra
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD25 Antibody [3C7] for use in flow cytometry, immunohistochemistry, and immunoprecipitation assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P01590
NCBI	P01590
Dilution Range	CD25 Antibody [3C7] for use in flow cytometry, immunohistochemistry, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



C57BL/6
mouse
cytotoxic T
lymphocyte
cel...