


Product Datasheet

Itgam Antibody (Biotin) (orb1251351)

Description	Itgam Antibody (Biotin)
Species/Host	Rat
Reactivity	Human, Mouse, Primate, Rabbit
Conjugation	Biotin
Tested Applications	ELISA, FC, IHC-Fr, NeA
Immunogen	C57BL/10 mouse splenic T cells and concanavalin A-activated C57BL/10 splenocytes
Target	Itgam
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD11b Antibody [M1/70] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and immunoprecipitation assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P05555
NCBI	P05555
Dilution Range	CD11b Antibody [M1/70] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.



BALB/c mouse splenocytes were stained wi...