

Product Datasheet

Itgal Antibody (Biotin) (orb1251359)

Description	Itgal Antibody (Biotin)
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC-Fr, NeA
Immunogen	Whole BW5147 cells
Target	Itgal
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD11a Antibody [I21/7] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and immunoprecipitation assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P24063
NCBI	P24063
Dilution Range	CD11a Antibody [I21/7] for use in flow cytometry, immunohistochemistry / immunocytochemistry, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.