

Product Datasheet

Itga4 Antibody (FITC) (orb1251397)

Description	Itga4 Antibody (FITC)
Species/Host	Rat
Reactivity	Human, Mouse
Conjugation	FITC
Tested Applications	ELISA, WB
Immunogen	P815 cells
Target	Itga4
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD49d Antibody [PS/2] for use in flow cytometry and immunohistochemistry assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q00651
NCBI	Q00651
Dilution Range	CD49d Antibody [PS/2] for use in flow cytometry and immunohistochemistry assays.
Expiration Date	12 months from date of receipt.