


## Product Datasheet

### Itga4 Antibody (orb1251398)

<b>Description</b>	Itga4 Antibody
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	P815 cells
<b>Target</b>	Itga4
<b>Preservatives</b>	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.5 mg/mL
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	CD49d Antibody [PS/2] for use in flow cytometry, immunohistochemistry, immunoprecipitation, blocking, and adhesion assays.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q00651</a>
<b>NCBI</b>	<a href="#">Q00651</a>
<b>Dilution Range</b>	CD49d Antibody [PS/2] for use in flow cytometry, immunohistochemistry, immunoprecipitation, blocking, and adhesion assays.
<b>Expiration Date</b>	12 months from date of receipt.



BALB/c mouse splenocytes were stained wi...