


Product Datasheet


Type IV Collagen Antibody (orb1251418)

Description	Type IV Collagen Antibody
Species/Host	Goat
Reactivity	Canine, Equine, Gallus, Guinea pig, Human, Mammal, Mouse, Porcine, Primate, Rabbit, Rat, Rodent, Sheep
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Type IV collagen
Preservatives	Supplied in borate buffered saline, pH 8.2; No preservatives or amine-containing buffer salts added
Form/Appearance	Liquid
Concentration	0.4 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Type IV Collagen Antibody for use in ELISA, immunohistochemistry / immunocytochemistry, electron microscopy, flow cytometry, western blotting, and immunoprecipitation assays.
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Type IV Collagen Antibody for use in ELISA, immunohistochemistry / immunocytochemistry, electron microscopy, flow cytometry, western blotting, and immunoprecipitation assays.
Expiration Date	12 months from date of receipt.




Human gastric cancer tissue was stained

...



Frozen sections of E12.5 eyes from Col4a...



Itga3-/-/Gpr56-/- E11.5 mouse neocortex ...