


Product Datasheet

IL17a Antibody (FITC) (orb1251459)

Description	IL17a Antibody (FITC)
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC
Immunogen	E. coli-expressed mouse IL-17
Target	IL17a
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	IL-17A Antibody [TC11-8H4] for use in flow cytometry assays.
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q62386
NCBI	Q62386
Dilution Range	IL-17A Antibody [TC11-8H4] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



PMA and Ionomycin-stimulated C57BL/6N mo...