


Product Datasheet

KRT19 Antibody (FITC) (orb1251463)

Description	KRT19 Antibody (FITC)
Species/Host	Mouse
Reactivity	Human
Conjugation	FITC
Tested Applications	Cell Separation, FC, IHC-Fr, IP
Immunogen	Recombinant C-terminal cytokeratin 19
Target	KRT19
Preservatives	Supplied in PBS/NaN3
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	Cytokeratin 19 Antibody [SB39g] for use in flow cytometry, immunocytochemistry, and western blotting assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P08727
NCBI	P08727
Dilution Range	Cytokeratin 19 Antibody [SB39g] for use in flow cytometry, immunocytochemistry, and western blotting assays.
Expiration Date	12 months from date of receipt.



Human pancreatic carcinoma cell line MIA...