


## Product Datasheet

### Mboat4 antibody (orb125149)

<b>Description</b>	Goat Polyclonal to Mboat4. NP_001119786.1, PMID: 19628676 and GenelD: 619373.
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Bovine, Human, Mouse, Porcine
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Target</b>	Mboat4 (mouse)
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	WB: Approx 50kDa band observed in Mouse Pancreas lysates (calculated MW of 49.9kDa according to NP_001119786.1). Recommended concentration: 0.5-2µg/ml. An additional band of unknown identity was also consistently observed at 38kDa. This band was successfully
<b>Protein Sequence</b>	CHLGLHYTEYYLGEP
<b>Clonality</b>	Polyclonal
<b>MW</b>	49.9
<b>Entrez</b>	<b>619373</b>
<b>Dilution Range</b>	ELISA: 1:16000, WB: 0.5-2 µg/ml
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of Mouse Pancreas ...