


Product Datasheet

CD8B Antibody (PE) (orb1251551)

Description	CD8B Antibody (PE)
Species/Host	Mouse
Reactivity	Gallus
Conjugation	PE
Tested Applications	Activity Assay, Blocking, Cytotoxic Assay, FC, IP
Immunogen	Chicken splenocytes
Target	CD8B
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD8beta Antibody [EP42] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	CD8beta Antibody [EP42] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



Chicken peripheral blood mononuclear cel...