

Product Datasheet

S-Tag Antibody [SBSTAGa] (PE) (orb1251576)

Description	S-Tag Antibody [SBSTAGa] (PE)
Species/Host	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	Recombinant S-tag protein
Preservatives	Supplied in PBS/NaN ₃ and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	S-Tag Antibody [SBSTAGa] for use in fluorescent assays.
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	S-Tag Antibody [SBSTAGa] for use in fluorescent assays.
Expiration Date	12 months from date of receipt.