


Product Datasheet

CD4 Antibody (PE/Cy5) (orb1251583)

Description	CD4 Antibody (PE/Cy5)
Species/Host	Mouse
Reactivity	Porcine
Conjugation	PE/Cy5
Tested Applications	Cell Separation, FC, IHC-Fr, IP
Immunogen	Fresh dd miniature swine thymocytes
Target	CD4
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD4 Antibody [74-12-4] for use in flow cytometry assays.
Isotype	IgG2b
Clonality	Monoclonal
Dilution Range	CD4 Antibody [74-12-4] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



Porcine peripheral blood lymphocytes wer...