

Product Datasheet

Icam2 Antibody (PE) (orb1251593)

Description	Icam2 Antibody (PE)
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE
Tested Applications	ELISA, FC, IHC-Fr, IP
Immunogen	COS cells transfected with mouse ICAM-2 cDNA and BW5147 AKR/J mouse thymoma cell line
Target	Icam2
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD102 Antibody [3C4] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P35330
NCBI	P35330
Dilution Range	CD102 Antibody [3C4] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.