


Product Datasheet

H2-K1 Antibody (PE) (orb1251596)

Description	H2-K1 Antibody (PE)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	PE
Tested Applications	FC
Immunogen	A.AL mouse splenocytes
Target	H2-K1
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	H-2Kk Antibody [36-7-5] for use in flow cytometry and immunohistochemistry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P01901
NCBI	P01901
Dilution Range	H-2Kk Antibody [36-7-5] for use in flow cytometry and immunohistochemistry assays.
Expiration Date	12 months from date of receipt.



AKR mouse splenocytes were stained with ...