


Product Datasheet

Cd19 Antibody (APC/Cy7) (orb1251761)

Description	Cd19 Antibody (APC/Cy7)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	APC/Cy7
Tested Applications	FC, ICC, IP, WB
Immunogen	CD19+ mouse pre-B cell line 300.19
Target	Cd19
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD19 Antibody [MB19-1] for use in flow cytometry assays.
Isotype	IgA
Clonality	Monoclonal
Uniprot ID	P25918
NCBI	P25918
Dilution Range	CD19 Antibody [MB19-1] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



BALB/c mouse splenocytes were stained wi...