


## Product Datasheet

### Cd19 Antibody (APC) (orb1251764)

<b>Description</b>	Cd19 Antibody (APC)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	APC
<b>Tested Applications</b>	ELISA, FC, ICC, IHC-P, IP, WB
<b>Immunogen</b>	CD19+ mouse pre-B cell line 300.19
<b>Target</b>	Cd19
<b>Preservatives</b>	Supplied in PBS/NaN3 and a stabilizing agent
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.1 mg/mL
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	CD19 Antibody [MB19-1] for use in flow cytometry assays.
<b>Isotype</b>	IgA
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P25918</b>
<b>NCBI</b>	<b>P25918</b>
<b>Dilution Range</b>	CD19 Antibody [MB19-1] for use in flow cytometry assays.
<b>Expiration Date</b>	12 months from date of receipt.



BALB/c mouse splenocytes were stained wi...