


Product Datasheet

Itgam Antibody (PE/Cy5) (orb1251776)

Description	Itgam Antibody (PE/Cy5)
Species/Host	Rat
Reactivity	Human, Mouse, Primate, Rabbit
Conjugation	PE/Cy5
Tested Applications	ELISA, FC, IHC-Fr, IP, WB
Immunogen	C57BL/10 mouse splenic T cells and concanavalin A-activated C57BL/10 splenocytes
Target	Itgam
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD11b Antibody [M1/70] for use in flow cytometry assays.
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P05555
NCBI	P05555
Dilution Range	CD11b Antibody [M1/70] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



BALB/c mouse splenocytes were stained wi...