


Product Datasheet

Itgam Antibody (PE/Cy7) (orb1251780)

Description	Itgam Antibody (PE/Cy7)
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE/Cy7
Tested Applications	ELISA, FC, ICC, IHC-Fr, NeA
Immunogen	Peritoneal macrophages from B6D2 hybrid mice
Target	Itgam
Preservatives	Supplied in PBS/NaN3 and a stabilizing agent
Form/Appearance	Liquid
Concentration	0.1 mg/mL
Storage	Store vial at 2-8°C
Note	For research use only
Application notes	CD11b Antibody [3A33] for use in flow cytometry assays.
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P05555
NCBI	P05555
Dilution Range	CD11b Antibody [3A33] for use in flow cytometry assays.
Expiration Date	12 months from date of receipt.



BALB/c mouse splenocytes were stained wi...