


## Product Datasheet

### Cd4 Antibody (APC/Cy7) (orb1251804)

<b>Description</b>	Cd4 Antibody (APC/Cy7)
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	APC/Cy7
<b>Tested Applications</b>	ELISA, FC, IHC-Fr, IHC-P, IP, Multiplex Assay, WB
<b>Immunogen</b>	Mouse CTL clone V4
<b>Target</b>	Cd4
<b>Preservatives</b>	Supplied in PBS/NaN3 and a stabilizing agent
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	0.1 mg/mL
<b>Storage</b>	Store vial at 2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	CD4 Antibody [GK1.5] for use in flow cytometry assays.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<b>P06332</b>
<b>NCBI</b>	<b>P06332</b>
<b>Dilution Range</b>	CD4 Antibody [GK1.5] for use in flow cytometry assays.
<b>Expiration Date</b>	12 months from date of receipt.



BALB/c mouse splenocytes were stained wi...