



Product Datasheet

CD31 antibody (orb12520)


Description	Rabbit polyclonal antibody to CD31
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human CD31/PECAM.
Preservatives	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	130kDa
Uniprot ID	P16284
Entrez	5175
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Rat placenta...



Immunohistochemical staining of Mouse placenta...



Immunohistochemical staining of Human placenta...