


Product Datasheet

TGM4 antibody (orb125201)

Description	Goat Polyclonal to TGM4. NP_003232.2, PMID: 20596668 and GeneID: 7047.
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	transglutaminase 4 (aa49-63)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	WB: Approx 75kDa band observed in lysates of cell line HeLa (calculated MW of 77.1kDa according to NP_003232.2). Recommended concentration: 1-3µg/ml.
Protein Sequence	NQPLQSYHQLKLEFS
Clonality	Polyclonal
MW	77.1
Entrez	7047
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa lysate using...