


Product Datasheet

HPV-16 Antibody [HPV16/1296] (orb1252157)

Description	HPV-16 Antibody [HPV16/1296]
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF, IHC-P
Immunogen	Human papilloma virus type 16 was used as the immunogen for the HPV-16 antibody.
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Immunofluorescence: 0.5-1 ug/mlIHC (FFPE): 0.5-1 ug/ml for 30 min at RTOptimal dilution of the HPV-16 antibody should be determined by the researcher.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Uniprot ID	NotKnown
Dilution Range	Immunofluorescence: 0.5-1 ug/mlIHC (FFPE): 0.5-1 ug/ml for 30 min at RTOptimal dilution of the HPV-16 antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.



IHC testing of FFPE human cervix with HPV...