



Product Datasheet

ACPP Antibody (orb1252197)

Description	ACPP Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC, WB
Immunogen	Human recombinant full length protein was used as the immunogen for this PSAP antibody.
Target	ACPP
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml IHC (FFPE): 0.5-1 ug/ml for 30 minutes at RT titering of the PSAP antibody may be required for optimal performance.
Isotype	IgG1, kappa
Clonality	Monoclonal
Uniprot ID	P15309
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml IHC (FFPE): 0.5-1 ug/ml for 30 minutes at RT titering of the PSAP antibody may be required for optimal performance.
Expiration Date	12 months from date of receipt.



IHC testing of human prostate carcinoma ...



SDS-PAGE Analysis of Purified, BSA-Free ...