

Product Datasheet

ALPI Antibody (orb1252201)

Description

ALPI Antibody

Species/Host

Mouse

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P

Applications

Immunogen

Recombinant human ALPL protein was used as the immunogen for the tissue-nonspecific Alkaline Phosphatase antibody.

Target

ALPI

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Immunofluorescence: 0.5-1 ug/ml
Immunohistochemistry (FFPE): 1-2 ug/ml for 30 min at RT
Prediluted format: incubate for 30 min at RT (1)
Optimal dilution of the Alkaline Phosphatase antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Isotype

IgG1, kappa

Clonality

Monoclonal

Uniprot ID


P09923

Dilution Range


Immunofluorescence: 0.5-1 ug/ml
Immunohistochemistry (FFPE): 1-2 ug/ml for 30 min at RT
Prediluted format: incubate for 30 min at RT (1)
Optimal dilution of the Alkaline Phosphatase antibody should be determined by the researcher.1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Expiration Date

12 months from date of receipt.



IHC testing of
FFPE ovarian
carcinoma
an...



IHC testing of
FFPE
testicular
carcinoma...