



Product Datasheet

Calnexin Antibody [CANX/1541] (orb1252359)


Description	Calnexin Antibody [CANX/1541]
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	A partial human recombinant protein was used as the immunogen for the Calnexin antibody.
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	WB: 0.5-1 ug/ml Flow Cytometry: 0.5-1ug/10 ⁶ cells Immunofluorescence: 0.5-1 ug/ml Immunohistochemistry (FFPE): 1-2 ug/ml for 30 min at RT Titration of the Calnexin antibody may be required for optimal performance.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Dilution Range	WB: 0.5-1 ug/ml Flow Cytometry: 0.5-1ug/10 ⁶ cells Immunofluorescence: 0.5-1 ug/ml Immunohistochemistry (FFPE): 1-2 ug/ml for 30 min at RT Titration of the Calnexin antibody may be required for optimal performance.
Expiration Date	12 months from date of receipt.



IHC testing of FFPE human small intestine...



IHC testing of FFPE human renal cell car...



IHC testing of FFPE human renal cell car...