

Product Datasheet

SLC3A2 Antibody (orb1252467)

Description	SLC3A2 Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF
Immunogen	Native leukocytes were used as the immunogen for the CD98 antibody.
Target	SLC3A2
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.2 mg/mL
Storage	Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml Optimal dilution of the CD98 antibody should be determined by the researcher.
Isotype	IgM
Clonality	Monoclonal
Uniprot ID	P08195
Dilution Range	Flow Cytometry: 0.5-1 ug/million cells in 0.1ml Immunofluorescence: 0.5-1 ug/ml Optimal dilution of the CD98 antibody should be determined by the researcher.
Expiration Date	12 months from date of receipt.