
Product Datasheet

HSPG2 Antibody (orb1252524)

Description

HSPG2 Antibody

Species/Host

Rat

Reactivity

Bovine, Fish, Human, Monkey, Mouse, Porcine

Conjugation

Unconjugated

Tested Applications

FC, IF, IHC-P, WB

Immunogen

A murine EHS laminin preparation was used as the immunogen for the Perlecan antibody.

Target

HSPG2

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
 Immunofluorescence: 0.5-1 ug/ml
 Immunohistochemistry (FFPE): 1-2 ug/ml for 30 min at RT
 (1) Prediluted format : incubate for 30 min at RT
 (2) Optimal dilution of the Perlecan antibody should be determined by the researcher.
 1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min
 2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Isotype


IgG2a, kappa

Clonality


Monoclonal

Uniprot ID
P98160
Dilution Range

Flow Cytometry: 0.5-1 ug/million cells in 0.1ml
 Immunofluorescence: 0.5-1 ug/ml
 Immunohistochemistry (FFPE): 1-2 ug/ml for 30 min at RT
 (1) Prediluted format : incubate for 30 min at RT
 (2) Optimal dilution of the Perlecan antibody should be determined by the researcher.
 1. Staining of formalin-fixed tissues requires boiling tissue sections in 10mM Citrate buffer, pH 6.0, for 10-20 min followed by cooling at RT for 20 min
 2. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section



Formalin-fixed, paraffin-embedded human ...



SDS-PAGE Analysis of Purified, BSA-Free ...