

Product Datasheet

EBV Early Antigens Antibody [1108-1] (orb1252664)

Description

EBV Early Antigens Antibody [1108-1]

Species/Host

Mouse

Reactivity

Virus

Conjugation

Unconjugated

Tested

IF

Applications
Immunogen

Affinity purified early antigen polypeptides from induced Raji cells precipitated with African Burkitt's lymphoma sera were used as the immunogen for this Early Antigens antibody.

Preservatives

PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide

Form/Appearance

Liquid

Concentration

0.2 mg/mL

Storage

Aliquot and Store at 2-8°C. Avoid freeze-thaw cycles.

Note

For research use only

Application notes

IF: 0.5-1 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Early Antigens antibody to be titrated up or down for optimal performance.

Isotype

IgG1, kappa

Clonality

Monoclonal

Dilution Range

IF: 0.5-1 ug/mlThe concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the Early Antigens antibody to be titrated up or down for optimal performance.

Expiration Date

12 months from date of receipt.