


Product Datasheet

NKIRAS1 Antibody (orb1252856)

Description	NKIRAS1 Antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A recombinant protein fragment from the central region of NKIRAS1 was used as the immunogen for this antibody.
Target	NKIRAS1
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Aliquot and Store at -20°C. Avoid freez-thaw cycles.
Note	For research use only
Application notes	Western blot: 5-10 ug/ml Provided assay concentrations are suggestions only, NKIRAS1 antibody titration may be required for optimal results.
Isotype	IgG1, kappa
Clonality	Monoclonal
Uniprot ID	Q9NYS0
Dilution Range	Western blot: 5-10 ug/ml Provided assay concentrations are suggestions only, NKIRAS1 antibody titration may be required for optimal results.
Expiration Date	12 months from date of receipt.



Western blot testing of human samples wi...