



Product Datasheet

OTUD5 Antibody (orb1252861)

Description	OTUD5 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A recombinant protein fragment from the C-terminal region of human OTUD5 was used as the immunogen for this antibody.
Target	OTUD5
Preservatives	PBS with 0.1 mg/ml rAlbumin and 0.05% sodium azide
Form/Appearance	Liquid
Concentration	0.5 mg/mL
Storage	Aliquot and Store at -20°C. Avoid freeze-thaw cycles.
Note	For research use only
Application notes	Western blot: 2-5 ug/ml Provided assay concentrations are suggestions only, antibody titration may be required for optimal results.
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96G74
Dilution Range	Western blot: 2-5 ug/ml Provided assay concentrations are suggestions only, antibody titration may be required for optimal results.
Expiration Date	12 months from date of receipt.



Western blot testing of human samples wi...



Western blot testing of mouse samples wi...