


Product Datasheet

PECAM1 Antibody (orb1252899)

Description	PECAM1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	PECAM-1 (phospho Tyr713) antibody was raised against a peptide sequence around phosphorylation site of Tyrosine 713 (T-V-Y (p)-S-E) derived from Human PECAM-1.
Target	PECAM1
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	82, 150 kDa
Uniprot ID	P16284
NCBI	NP_000433.4
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from J...