


Product Datasheet

ABL1 Antibody (orb1252904)

Description	ABL1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	c-Abl (phospho Tyr245) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 245 (T-V-Y (p)-G-V) derived from Human c-Abl.
Target	ABL1
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	135 kDa
Uniprot ID	P00519
NCBI	NP_005148.2
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from C...