
Product Datasheet

CAMK4 Antibody (orb1252912)

Description

CAMK4 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

CaMK4 (phospho Thr196/200) antibody was raised against a peptide sequence around phosphorylation site of threonine 200 (M-K-T (p)-V-C) derived from Human CaMK4 or threonine 196 (M-K-T (p)-V-C) derived from Mouse CaMK4.

Target

CAMK4

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000

Isotype

IgG

Clonality

Polyclonal

MW


60 kDa

Uniprot ID
[Q16566](#)
NCBI
[NP_001735.1](#)
Dilution Range

Western Blot: 1:500~1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis of extracts from K...