


## Product Datasheet

### BLNK Antibody (orb1252929)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | BLNK Antibody   |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse  |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | WB  |
| <b>Immunogen</b>           | BLNK (phospho Tyr96) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 96 (D-S-Y (p)-E-P) derived from Human BLNK .  |
| <b>Target</b>              | BLNK  |
| <b>Preservatives</b>       | Antibody supplied in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Concentration</b>       | 1 mg/mL   |
| <b>Storage</b>             | Store antibody at -20°C for up to one year.   |
| <b>Note</b>                | For research use only   |
| <b>Application notes</b>   | Western Blot: 1:500~1:1000  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 50 kDa  |
| <b>Uniprot ID</b>          | <a href="#">Q8WV28</a>  |
| <b>NCBI</b>                | <a href="#">NM_001114094.1</a>  |
| <b>Dilution Range</b>      | Western Blot: 1:500~1:1000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Western blot analysis of extracts from C...