


Product Datasheet

VIM Antibody (orb1252952)

Description	VIM Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Vimentin (phospho Ser56) antibody was raised against a peptide sequence around phosphorylation site of serine 56 (A-S-S (p)-P-G) derived from Human Vimentin.
Target	VIM
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	P08670
NCBI	NP_003371.2
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from A...