


Product Datasheet

IL10RA Antibody (orb1252968)

Description	IL10RA Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	IL-10R alpha (phospho Tyr496) antibody was raised against a peptide sequence around phosphorylation site of Tyrosine496 (K-G-Y (p)-L-K) derived from Human IL-10R alpha.
Target	IL10RA
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	63 kDa
Uniprot ID	Q13651
NCBI	NP_001549.2
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts HUVEC ...