


Product Datasheet

NFKB2 Antibody (orb1252971)

Description	NFKB2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	NF-κB p100 (phospho Ser872) antibody was raised against a peptide sequence around phosphorylation site of serine872 (S-Q-S (p)-V-E) derived from Human NF-κB p100 .
Target	NFKB2
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	100 kDa
Uniprot ID	Q00653
NCBI	NP_001070962.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from R...