


Product Datasheet

MTOR Antibody (orb1252975)

Description	MTOR Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	mTOR (phospho Thr2446) antibody was raised against a peptide sequence around phosphorylation site of threonine2446 (T-R-T (p)-D-S) derived from Human mTOR .
Target	MTOR
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	220 kDa
Uniprot ID	P42345
NCBI	NP_004949.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from N...