


Product Datasheet

NTRK1 Antibody (orb1252986)

Description	NTRK1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Trk A (phospho Tyr680/681) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 680/tyrosine 681 (T-D-Y (p)-Y (p)-R-V) derived from Human Trk A.
Target	NTRK1
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	87 kDa
Uniprot ID	P04629
NCBI	NP_001007793.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from Jurkat starved cells.