


Product Datasheet

MAP3K7 Antibody (orb1252990)

Description	MAP3K7 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	MAP3K7 (phospho Thr187) antibody was raised against a peptide sequence around phosphorylation site of threonine 187 (H-M-T (p)-N-N) derived from Human MAP3K7.
Target	MAP3K7
Preservatives	Antibody supplied in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Form/Appearance	Liquid
Concentration	1 mg/mL
Storage	Store antibody at -20°C for up to one year.
Note	For research use only
Application notes	Western Blot: 1:500~1:1000
Isotype	IgG
Clonality	Polyclonal
MW	60 kDa
Uniprot ID	O43318
NCBI	NP_003179.1
Dilution Range	Western Blot: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from N...