
Product Datasheet

TUBA1A, TUBA1B, TUBA1C Antibody (orb1253004)

Description

TUBA1A, TUBA1B, TUBA1C Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

TUBA1/3/4 (phospho Tyr272) antibody was raised against a peptide sequence around phosphorylation site of tyrosine 272 (A-T-Y (p)-A-P) derived from Human TUBA1/3/4.

Target

TUBA1A, TUBA1B, TUBA1C

Preservatives

Antibody supplied in phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Form/Appearance

Liquid

Concentration

1 mg/mL

Storage

Store antibody at -20°C for up to one year.

Note

For research use only

Application notes

Western Blot: 1:500~1:1000

Isotype

IgG

Clonality

Polyclonal

MW

50-55 kDa

Uniprot ID

[Q6PEY2](#), [Q13748](#), [Q9BQE3](#), [P68363](#), [Q71U36](#), [P68366](#)

NCBI

[NP_001257328.1](#)

Dilution Range

Western Blot: 1:500~1:1000

Expiration Date

12 months from date of receipt.

